


Kiato



Long Stitch Cutter • Standard Stitch Cutter • Mini Stitch Cutter

After more than 3 decades of supremacy in the Indian market, the Kiato range of surgical blades, disposable scalpels and stitch cutters made their global presence felt in 2000. Today, over 80 countries bear testimony to the excellent quality and superior performance of our products, both under the Kiato brand name as well as in private label.





Long Stitch Cutter



Ref 11603



Ref 11605

Individually packed in blister pack. Also can be supplied in bulk packs for kits. The long stitch cutter offers ease of use with better visibility and safe handling.

Material: Swedish stainless steel 100 per box
Gamma Sterile: 2.5 M Rad

Standard Stitch Cutter



Ref 11600 Individually packed in easy peel foil with VCI paper protection. This is traditionally the most commonly preferred stitch cutter competitively priced.

Material: Swiss high carbon steel 100 per box
Gamma Sterile: 2.5 M Rad

Mini Thumb Stitch Cutter



Mini Stitch Cutter



Ref 11602 Individually packed in easy peel foil. Also offered in bulk packing assembled on a thumb scalpel with a safety cap. This offers a convenient alternative to the conventional long stitch cutter.

Material: Swedish stainless steel 100 per box
Gamma Sterile: 2.5 M Rad

All stitch cutters have a rounded tip which makes them ideal for safely cutting sutures and minimizing patient discomfort. We offer a full range competitively priced for better value.

Products are not shown to scale



All products are CE marked as per the European Medical Device Directive 93/42/EEC.



MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India

Tel. +91 124 47 64 900
Fax +91 124 47 64 911

info@adityadisposed.com
www.adityadisposed.com

A subsidiary of Kehr Surgical Pvt. Ltd.